

Baby Bath Promotional Thermometer

SAFE
ACCURATE
UNBREAKABLE
EASY TO USE
DURABLE

With custom printing that conveys your message as well as heightening the understanding of baby bath temperature, this thermometer is the ideal promotional gift to mothers and expecting mothers. It is easy to include in a mailing or a gift basket.

The thermometer with your message is long lived, because even when they get older the children will insist on using it every time they have a bath.

Reading:

The liquid crystals light up when their respective temperatures are reached. Read the brightest number. If two numbers are equally bright, the temperature is in between.


Instructions:

Stir water to get even temperature.

Put thermometer fully into water and keep it there for 15 seconds. Temperature must be read while thermometer is still under water.

Care of thermometer:

Do not leave in direct sunshine. Wipe dry after use.



Specifications:

Temperature range	+30 +50°C	Dimensions	146 x 66mm
Resolution	2°C increments	Weight	3.8 grams
Accuracy	± 0.5°C		

Temperature Technology

263 Gilbert Street

Adelaide SA 5000

Tel: 08-8231 1266

tec.com.au

Fax 08-8231 1212

email: sales@t-

Web: www.t-tec.com.au